

Sharpsafe 0.2 & 0.3 Litre



0.2 Litre

0.3 Litre

Features:

- > Welded construction for added strength & safety
- > Lid opening with simple pressure with just a single hand
- > Flip top lid for temporary and final closure
- > In-mould label which cannot be removed, even after autoclaving
- > Handle for a safe transport
- > Puncture resistant base
- > Autoclavable to 135°C for 20 minutes for an empty and open container
- > Conforms to Australian, British and European standards

Specifications:

Size	0.2L	0.3L
Product code:	36-4173	36-4175
Material:	Polypropylene	Polypropylene
Dimensions – Container (mm):	79 x 37 x 146	105 x 50 x 130
Dimensions – Opening (mm):	25 x 20	40 x 30
Capacity to Fill Line:	170mL	255mL
Weight:	55g	64g
Packaging		
Quantity per box:	100	100
Dimensions:	398 x 321 x 306	445 x 274 x 506
Weight:	6.35kg	8kg

Sharpsafe 0.45 & 0.6 Litre



0.45 Litre



0.6 Litre

Features:

- > Fully assembled, welded construction for added strength & safety
- > Lid opening with simple pressure with just a single hand
- > Flip top lid for temporary & final closure
- > In-mould label which cannot be removed, even after autoclaving
- > Fold away handle for a safe transport
- > Puncture resistant base
- > Autoclavable to 135°C for 20 minutes for an empty & open container
- > Conforms to Australian, British & European standards

Specifications:

Size	0.45L	0.6L
Product code:	36-4170	36-4150
Material:	Polypropylene	Polypropylene
Dimensions – Container (mm):	103 x 50 x 169	106 x 54 x 221
Dimensions – Opening (mm):	40 x 30	40 x 30
Capacity to Fill Line:	360mL	500mL
Weight:	88g	110g
Packaging		
Quantity per box:	100	100
Dimensions:	503 x 443 x 342	530 x 455 x 460
Weight:	9.2kg	12kg

Sharpsafe 1 Litre



1 Litre

Features:

- > Safe, quick and simple 4 pressure assembly with 4 audible « click »
- > 2 openings: Universal removal & unobstructed opening
- > 2 closure positions on the flip top lid: temporary & final closure
- > Reopening of the temporary closure with a simple finger pressure on the lid
- > In-mould label which cannot be removed, even after autoclaving
- > Handle for safe transport
- > Autoclavable to 135°C for 20 minutes for an empty, open container
- > Conforms to Australian, British & European standards

Specifications:

Size

1L

Product code:	36-4160
Material:	Polypropylene
Dimensions – Container (mm):	108 x 90 x 180
Dimensions – Opening (mm):	55 x 35
Capacity to Fill Line:	850mL
Weight:	126g

Packaging

Quantity per box:	100
Dimensions:	545 x 390 x 430
Weight:	14kg

Sharpsafe 2 & 4 Litre



2 Litre

4 Litre

Features:

- > Safe, quick and simple 4 pressure assembly with 4 audible « click »
- > Wide opening & large headroom providing users with sufficient safety buffer
- > Disconnection for all types of needles, even for scalpel blades & insulin pens
- > Internal safety flap for extra protection
- > 2 closure positions on the flip top lid: temporary & final closure
- > Handle for safe transport
- > Accessories: Stands, hooks & brackets available
- > Autoclavable to 135°C for 20 minutes for an empty, open container
- > Conforms to Australian, British & European standards

Specifications:

Size	2L	4L
Product code:	36-4140	36-4100
Material:	Polypropylene	Polypropylene
Dimensions – Container (mm):	195 x 120 x 178	175 x 175 x 245
Dimensions – Opening (mm):	55 x 50	55 x 50
Capacity to Fill Line:	1.7L	3.2L
Weight:	220g	310g
Packaging		
Quantity per box:	50	50
Dimensions:	585 x 385 x 385	843 x 355 x 420
Weight:	12kg	17kg

Sharpsafe 6 & 7 Litre



Features:

- > Safe, quick and simple 4 pressure assembly with 4 audible « click »
- > Wide opening & large headroom providing users with sufficient safety buffer
- > Disconnection for all types of needle even for the scalpel blade & insulin pen
- > Internal safety flap for extra protection
- > 2 closure positions on the flip top lid: temporary & final closure
- > Handle for safe transport
- > Accessories: Stands, hooks & brackets available
- > Autoclavable to 135°C for 20 minutes for an empty, open container
- > Conforms to Australian, British & European standards

Specifications:

Size	6L	7L
Product code:	36-4206	36-4110
Material:	Polypropylene	Polypropylene
Dimensions – Container (mm):	175 x 175 x 300	180 x 180 x 328
Dimensions – Opening (mm):	100 x 50	105 x 50
Capacity to Fill Line:	4.8L	5.8L
Weight:	382g	350g
Packaging		
Quantity per box:	40	40
Dimensions:	880 x 365 x 400	800 x 365 x 415
Weight:	17kg	16kg

Sharpsafe 8/12 & 24 Litre



Features:

- > Safe, quick and simple 4 pressure assembly with 4 audible « click »
- > Wide opening & large headroom providing users with sufficient safety buffer
- > Disconnection for all types of needle even for the scalpel blade & insulin pen
- > Internal safety flap for extra protection
- > 2 closure positions on the flip top lid: temporary & final closure
- > Handle for safe transport
- > Autoclavable to 135°C for 20 minutes for an empty, open container
- > Conforms to British & European standards

Specifications:

Size	8L	12L	24L
Product code:	36-4208	36-4212	36-4224
Material:	Polypropylene	Polypropylene	Polypropylene
Dimensions – Container (mm):	289 x 189 x 240	289 x 189 x 310	450 x 230 x 335
Dimensions – Opening (mm):	170 x 55	170 x 55	365 x 95
Capacity to Fill Line:	6.5L	9.8L	18.7L
Weight:	520g	597g	1.1kg
Packaging			
Quantity per box:	20	20	15
Dimensions:	590 x 410 x 440	590 x 410 x 490	445 x 445 x 795
Weight:	12kg	12.5kg	18.5kg

Protective Access Range 9 & 13 Litre



9 Litre PA



13 Litre PA



Protected Access

Features:

- > Protected Access allows easy disposal of sharps but provides added security for both the user & the general public.
- > Protected Access insert prevents risk of needle stick injury
- > Safe, quick & simple 4 pressures assembly with 4 audible « click »
- > Extra wide opening with horizontal fill area to easily & safely dispose all sharps
- > Handle for a safe transport
- > In-mould label – offers permanent audit trail which cannot be removed
- > Autoclavable to 135°C for 20 minutes for an empty, open container
- > Conforms to Australian, British & European standards

Specifications:

Size	9L PA	13L PA
Product code:	36-4117	36-4115
Material:	Polypropylene	Polypropylene
Dimensions – Container (mm):	195 x 292 x 291	195 x 292 x 355
Dimensions – Opening (mm):	202 x 38	202 x 38
Capacity to Fill Line:	7.5L	10.4L
Weight:	619g	709g
Packaging		
Quantity per box:	20	20
Dimensions:	585 x 385 x 530	590 x 400 x 590
Weight:	13.8kg	15.5kg

Clinisafe Cardboard Cartons



Features:

The Clinisafe range of Cardboard Clinical Waste Containers is constructed from recycled cardboard and contains an integral polyethylene liner with ties for closure.

The environmentally friendly, lightweight design allows the Clinisafe Cardboard Clinical Waste Container range to be delivered folded flat prior to use reducing transport costs and minimising on storage space.

Ideal for solid and soft clinical waste.

This range is approved to UN3291 and NFX 30-507.

Specifications:

Size	12L	25L	50L
Product code:	36-3120	36-3250	36-3500
Material:	197 x 197 x 378	264 x 264 x 370	293 x 293 x 575
Dimensions – Container (mm):	200 x 200 x 385	267 x 267 x 380	295 x 295 x 580
Dimensions – Opening (mm):	60 microns	60 microns	60 microns
Capacity to Fill Line:	11L	23L	48L
Weight:	177g	298g	495g
Packaging			
Shipping/storage dimensions:	55 x 540 x 390	70 x 530 x 570	70 x 790 x 585
Quantity per batch:	10	10	10
Quantity per pallet:	800	800	400

Features & Benefits




74% Recycled Cardboard

12L is 65% recycled
25L is 70% recycled
50L is 74% recycled




Ultra Lightweight

Lighter than traditional plastic containers, lowering transport costs



Emissions are minimised

Transport costs reduced



Stacking

Stackable when full to convenient storage



FSC approved materials

from sustainable sources

Indicating the product has been manufactured using board from sustainable forestry sources



Bar Coding

Unique barcode for traceability

All products also contain a unique identifier bar code which can allow single unit traceability



Accessories

25L Plinth **44100004**
for use with all 25L cardboard containers

50L Short Plinth **44100001**
for use with all 50L short cardboard containers

50L Tall Plinth **44100005**
for use with 50L yellow cardboard container



Storage

Cardboard cartons are available as pallets of 400 (25L & 50L); pallets of 800 (12L); or as banded 10's.

Flatpack solution to save on crucial stock and transport space.

Products in banded 10's are easily stored and intended for use in areas of low usage.



Stacking

Cartons of the same size can be efficiently stacked, minimising pre-disposal storage space.

Clinisafe Fluid Tight Containers



Full Range



Stand to Closure

Features:

Ranging from 3 litres to 11 litres in capacity the Clinisafe Plastic Clinical Waste Containers have been designed to hold clinical, anatomical & laboratory waste.

- > Fluid leak resistant seal with temporary closure feature
- > In-mould label – offers permanent audit trail which cannot be removed.
- > Simple 4 Click Final Closure
- > Wall or trolley mountable for extra stability
- > Carry Handle
- > Interstackable for efficient space utilisation when lid is open or stackable when closed
- > Complies with UN 3291 and BS EN ISO 23907:2012 standard

Specifications:

Size	3L	5L	8L	11L
Product code:	36-3030	36-3050	36-3080	36-3110
Material:	Polypropylene	Polypropylene	Polypropylene	Polypropylene
Dimensions - Container:	174 x 174 x 185	174 x 174 x 322	315 x 230 x 250	315 x 230 x 320
Capacity to Fill Line:	2.22L	4.25L	6.1L	9L
Weight:	255g	388g	537g	622g
Packaging				
Quantity per box:	20	20	20	20
Dimensions:	353 x 340 x 338	330 x 345 x 585	465 x 330 x 735	630 x 314 x 783
Weight:	5.7kg	8.9kg	10.7kg	13.6kg

Sharpsafe Accessories

Brackets, Stands & Trays



Features:

- > Light weight, discrete and secure.
- > Universal wall clips
- > Universal loop bracket

Wall Loop Bracket



Features:

- > Stainless steel
- > Universal wall clips
- > Universal loop bracket


Stands & Trays



Features:


- > The Sharpsafe Tray forms an important part of the Sharpsafe Near Patient Sharps Disposal System (NPSD™) and has been designed to provide a large, easy to clean procedure tray.
- > The Sharpsafe Tray has been designed to accommodate the 1, 2 or 3L Sharpsafe containers
- > Stands for Sharpsafe containers enhance stability, improving safety, security & prevent accidental spillage or tipping

Assembly, Use & Closure



1. Assembly

Place the base of the container on a firm surface, then attach the lid over the base and firmly apply pressure to each corner so that it engages with a “click”.




2. Temporary Close

Ensure temporary closure is engaged when the container is not in use.



3. Re-Open

Push tab up to release lid from temporary closure.



4. Final Closure

Place the container onto suitable surface to lock the Sharpsafe container ready for disposal.

Complete Range Summary					
Size	Cat. No.	Unit Dimensions	Carton Dimensions	Unit Qty per Carton	Carton Weight
0.2L	36-4173	146 x 79 x 37	310 x 398 x 330	100	6.5kg
0.3L	36-4175	130 x 105 x 50	275 x 510 x 445	100	8kg
0.45L	36-4170	169 x 103 x 50	350 x 505 x 445	100	9kg
0.6L	36-4150	221 x 106 x 54	465 x 540 x 460	100	12kg
1L	36-4160	180 x 108 x 90	430 x 545 x 390	100	14kg
2L	36-4140	178 x 195 x 120	385 x 585 x 385	50	12kg
4L	36-4100	245 x 175 x 175	417 x 843 x 355	50	17kg
6L	36-4206	300 x 175 x 175	880 x 365 x 400	40	17kg
7L	36-4110	328 x 180 x 180	800 x 415 x 365	40	16kg
8L	36-4208	240 x 289 x 189	590 x 410 x 440	20	12kg
12L	36-4215	310 x 289 x 189	590 x 410 x 490	20	13kg
24L	36-4224	335 x 450 x 230	445 x 445 x 795	15	19kg
9L PA	36-4117	292 x 195 x 291	585 x 385 x 530	20	14kg
13L PA	36-4115	355 x 195 x 291	590 x 400 x 590	20	16kg
12L CC	36-3120	385 x 200 x 200	110 x 400 x 595	10	2kg
25L CC	36-3250	380 x 267 x 267	110 x 528 x 595	10	3kg
50L CC	36-3500	580 x 295 x 295	110 x 595 x 780	10	5kg
3L CF	36-3030	185 x 174 x 174	353 x 340 x 338	20	6kg
5L CF	36-3050	322 x 174 x 174	330 x 345 x 585	20	9kg
8L CF	36-3080	250 x 315 x 230	465 x 330 x 735	20	11kg
11L CF	36-3100	320 x 315 x 230	630 x 314 x 783	20	14kg
Accessories		Product Code		Qty/Pack	
4/6/7 Litre Wall Bracket		36-4340		5	
2/4/7 Litre Clip Bracket		36-4310		10	
0.3/0.45/0.6 Litre Loop Bracket		36-4320		10	
1/2 Litre Tray		36-4300		10	
1/2 Litre Stand		36-4150B		10	
9/13 Litre Stand		36-4117B		10	